

APB486Hu61 100µg

Active Calreticulin (CALR)

Organism Species: *Homo sapiens* (Human)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Eukaryotic expression.

Host: 293F cell

Residues: Glu18~Leu417

Tags: N-terminal His-tag

Purity: >98%

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).

Buffer Formulation: PBS, pH7.4, containing 5% Trehalose .

Original Concentration: 200µg/mL

Applications: Cell culture; Activity Assays.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 4.3

Predicted Molecular Mass: 48.1kDa

Accurate Molecular Mass: 60kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```
EPAVYFKEQFLDGDGWTSRWIESKHKSDFGKFLVSSGKFYGDEEKDKGLQTSQDARFYALSASFEPFS
NKGQTLVVQFTVKHEQNIIDCGGGYVKLFPNSLDQTDMHGDSEYNIMFGPDI CGPGTKKVHVIFNYKGG
NVLINKDIRCKDDEFTHLYTLIVRPDNTYEVKIDNSQVESGSLEDDWDFLPPKKIKDPDASKPEDWDE
RAKIDDPTDSKPEDWDKPEHIPDPDAKKPEDWDEEMDGEWEPPVIQNPEYKGEWKPRQIDNPDYKGTW
IHPEIDNPEYSPDPSIYAYDNFGVLGLDLWQVKSGTIFDNFLITNDEAYAEFEFGNETWGVTKAAEKQM
KDKQDEEQRLKEEEEKDKRKEEEEAEDKEDDEDKDEDEEDEDKEEDEEEDVPGQAKDEL
```

[ACTIVITY]

Calreticulin (CALR) is an endoplasmic reticulum (ER)-resident protein involved in a spectrum of cellular processes. In healthy cells, CALR operates as a chaperone and Ca²⁺ buffer to assist correct protein folding within the ER. Calreticulin is also found in the nucleus, suggesting that it may have a role in transcription regulation. Systemic lupus erythematosus is associated with increased autoantibody titers against calreticulin. Recurrent mutations in calreticulin have been linked to various neoplasms, including the myeloproliferative type. Adhesive glycoprotein THBS1 has been shown to play roles in platelet aggregation, angiogenesis and tumorigenesis, the protein is an interactor of CALR, thus a functional binding ELISA assay was conducted to detect the interaction of recombinant human CALR and recombinant rat THBS1. Briefly, biotin-linked CALR were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 ul were then transferred to THBS1-coated microtiter wells and incubated for 1h at 37°C. Wells were washed with PBST 3 times and incubation with Streptavidin-HRP for 30min, then wells were aspirated and washed 5 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C.

Finally, add 50 µl stop solution to the wells and read at 450 nm immediately. The binding activity of CALR and THBS1 was shown in Figure 1, the EC50 for this effect is 0.018 µg/mL.




Figure 1. The binding activity of recombinant human CALR and recombinant rat THBS1

[IDENTIFICATION]




Figure 1. Gene Sequencing (extract)



Figure 2. SDS-PAGE

Sample: Active recombinant CALR, Human

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.