

PAA665Hu71

Biotin-Linked Antibody to Cross Linked C-Telopeptide Of Type I Collagen (CTXI)

**Organism Species: Homo sapiens (Human)** 

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

## [ PRODUCT INFORMATION ]

Immunogen: CTXI-OVA Clonality: Polyclonal Conjugation: Biotin

Host: Rabbit

Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 200µg/mL

**UOM**: 5μg

## [ IMMUNOGEN INFORMATION ]

Immunogen: Synthetic Peptide, CTXI conjugated to OVA.

Accession No.: CPA665Hu21

**Sequence:** The target peptide sequence is listed below.

**EKAHDGGRYYRA** 

## [ANTIBODY SPECIFITY]

The antibody is a rabbit polyclonal antibody raised against CTXI. It has been selected for its ability to recognize CTXI in immunohistochemical staining and western blotting.



## [APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200 Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

# [CONTENTS]

**Form & Buffer:** Supplied as solution form in PBS, pH7.4, containing 0.02% NaN<sub>3</sub>, 50% glycerol.

# [STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.