

RPF052Hu01 100ug Recombinant Tripartite Motif Containing Protein 5 (TRIM5) Organism Species: Homo sapiens (Human) *Instruction manual* 

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

## [PROPERTIES]

Residues: Ala2~Thr261 **Tags:** Two N-terminal Tags, His-tag and T7-tag kDa 70 Accession: Q9C035 44 Host: E. coli 33 Subcellular Location: Cytoplasm. 26 **Purity:** >90% 22 Endotoxin Level: <1.0EU per 1µg 18 (determined by the LAL method). Formulation: Supplied as lyophilized form in 20mM Tris, 14 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 10 0.01% sarcosyl, 5% trehalose, and preservative. 15% SDS-PAGE Predicted isoelectric point: 5 Predicted Molecular Mass: 33.7kDa Applications: SDS-PAGE; WB; ELISA; IP. (May be suitable for use in other assays to be determined by the end user.)

## [ <u>USAGE</u> ]

Reconstitute in sterile ddH<sub>2</sub>O.



## [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## [<u>SEQUENCES</u>]

The sequence of the target protein is listed below.

ASGILVNVK EEVTCPICLE LLTQPLSLDC GHSFCQACLT ANHKKSMLDK GESSCPVCRI SYQPENIRPN RHVANIVEKL REVKLSPEGQ KVDHCARHGE KLLLFCQEDG KVICWLCERS QEHRGHHTFL TEEVAREYQV KLQAALEMLR QKQQEAEELE ADIREEKASW KTQIQYDKTN VLADFEQLRD ILDWEESNEL QNLEKEEEDI LKSLTNSETE MVQQTQSLRE LISDLEHRLQ GSVMELLQGV DGVIKRTENV T