

SPA957Hu02 50µg

Synthetic Procollagen I N-Terminal Propeptide (PINP)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Synthetic Peptide

Residues: Thr78-Pro93

Purity: >95%

Traits: Freeze-dried powder

Buffer Formulation: PBS, pH7.4.

Original Concentration: 5mg/mL

Predicted isoelectric point: 4.3

Predicted Molecular Mass: 1633.8Da

Applications: Immunogen; Blocking Peptide.

[USAGE]

Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

TKNCPGAEVP EGECCP

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.