

APD862Hu01 100µg

Active Tyrosine Aminotransferase (TAT)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Apr. 2016)

[PROPERTIES]

Source: Prokaryotic expression.

Host: E. coli

Residues: Cys221~Lys454
Tags: N-terminal His-tag

Purity: >96%

Buffer Formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.05% sarcosyl

and 5% trehalose.

Applications: Cell culture; Activity Assays.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.5

Predicted Molecular Mass: 30.0kDa

Accurate Molecular Mass: 33kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE] CGSVFSKRHL QKILAVAARQ CVPILADEIY GDMVFSDCKY EPLATLSTDV PILSCGGLAK RWLVPGWRLG WILIHDRRDI FGNEIRDGLV KLSQRILGPC TIVQGALKSI LCRTPGEFYH NTLSFLKSNA DLCYGALAAI PGLRPVRPSG AMYLMVGIEM EHFPEFENDV EFTERLVAEQ SVHCLPATCF EYPNFIRVVI TVPEVMMLEA CSRIQEFCEQ HYHCAEGSQE ECDK

[ACTIVITY]

Tyrosine aminotransferase (TAT) is an enzyme present in the liver and catalyzes the conversion of tyrosine to 4-hydroxyphenylpyruvate. In humans, the tyrosine aminotransferase protein is encoded by the TAT gene. A deficiency of the enzyme in humans can result in what is known as Type II Tyrosinemia, wherein there is an abundance of tyrosine as a result of tyrosine failing to undergo an aminotransferase reaction to form 4-hydroxyphenylpyruvate. Aminotransferase as a dimer has two identical active sights. Lys280 is attached to PLP, which is held in place via two nonpolar amino acid side chains; phenylalanine and isoleucine (see thumbnail on right). The PLP is also held in place by hydrogen bonding to surrounding molecules mainly by its phosphate group. Besides, Heat Shock 70kDa Protein 8 (HSPA8) has been identified as an interactor of TAT, thus a binding ELISA assay was conducted to detect the interaction of recombinant human TAT and recombinant human HSPA8. Briefly, TAT were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100uL were then transferred to HSPA8-coated microtiter wells and incubated for 2h at 37 °C. Wells were washed with PBST and incubated for 1h with anti-TAT pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated

Cloud-Clone Corp.

and washed 3 times. With the addition of substrate solution , wells were incubated 15-25 minutes at $37\,^{\circ}$ C. Finally, add 50μ L stop solution to the wells and read at 450nm immediately. The binding activity of TAT and HSPA8 was shown in Figure 1, and this effect was in a dose dependent manner.




Figure 1. The binding activity of TAT with HSPA8.

[IDENTIFICATION]




Figure 2. SDS-PAGE

Sample: Active recombinant TAT, Human

Cloud-Clone Corp.




Figure 3. Western Blot

Sample: Recombinant TAT, Human;

Antibody: Rabbit Anti-Human TAT Ab (PAD862Hu01)