
#### PAA726Bo01 Polyclonal Antibody to Transthyretin (TTR) Organism Species: Bos taurus; Bovine (Cattle) *Instruction manual*

#### FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES



#### [IMMUNOGEN INFORMATION]

Immunogen: Recombinant TTR (Ser22~Ala147) expressed in *E.coli*.

#### USCN Accession No.: RPA726Bo01

**Sequence:** The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHSGSEF-SVGAGEPKC PLMVKVLDAV RGSPAANVGV KVFKKAADET WEPFASGKTS ESGELHGLTT EDKFVEGLYK VELDTKSYWK SLGISPFHEF AEVVFTANDS GPRHYTIAAL LSPYSYSTTA LVSSPKA

## [ANTIBODY SPECIFITY]

The antibody is a rabbit polyclonal antibody raised against TTR. It has been selected for its ability to recognize TTR in immunohistochemical staining and western blotting.

### [APPLICATIONS]

Western blotting: 1:100-400 Immunocytochemistry in formalin fixed cells: 1:100-500 Immunohistochemistry in formalin fixed frozen section: 1:100-500 Immunohistochemistry in paraffin section: 1:50-200 Enzyme-linked Immunosorbent Assay: 1:100-200 Optimal working dilutions must be determined by end user.

# [<u>CONTENTS</u>]

**Form & Buffer:** Supplied as solution form in PBS, pH7.4, containing 0.02% NaN<sub>3</sub>, 50% glycerol.

## [ QUALITY CONTROL ]

**Content:** The quality control contains recombinant TTR (Ser22~Ala147) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemilumescent (ECL). **Note:** The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading Buffer: 100mM Tris(pH8.8), 2% SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, NaN $_3$  0.02%.

# [<u>STORAGE</u>]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.