




P91343Po01 Bcl2 Associated X Protein (Bax) Organism: Sus scrofa; Porcine (Pig) *Instruction manual* 

# FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

## [ DESCRIPTION ]



Unique product Superb quality Client favorite Nicest service Q ISO9001:2008; Q ISO13485:2003;  $\zeta \in$ 





## [PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

#### [STORAGE AND STABILITY]

#### Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

### [SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below. MGHHHHHHSGSEF- DGSGEQPRG GGPTSSEQIM KTGALLLQGF IQDRAGRMGG ETPELGLEQV PQDASTKKLS ECLKRIGDEL DSNMELQRMI AAVDTDSPRE VFFRVAAEMF ADGNFNWGRV VALFYFASKL VLKALCTRVP ELIRTIMGWT LDFLRDRLLG WIQDQGGWDG LLSYFGTPTW QTVTIFVAGV LTASLTIWKK MG

Unique product Superb quality Client favorite Nicest service @ ISO9001:2008; @ ISO13485:2003;  $\in$