

P93079Mu01
Carbonic Anhydrase XII (CA12)
Organism: Mus musculus (Mouse)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
 NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)


[DESCRIPTION]

Protein Names: Carbonic Anhydrase XII
Synonyms: CA12, Car12
Species: Mouse
Size: 100µg
Source: *Escherichia coli*-derived
Subcellular Location: Membrane; Single-pass type I membrane protein.

[PROPERTIES]

Residues: Ala161~Ile318 (Accession # Q8CI85), with N-terminal His-Tag.
Grade & Purity: >95%, 21kDa as determined by SDS-PAGE reducing conditions.
Formulation: Supplied as lyophilized form in PBS, pH 7.4, containing 5% sucrose, 0.01% sarcosyl.
Endotoxin Level: <1.0 EU per 1µg (determined by the LAL method).
Applications: SDS-PAGE; WB; ELISA; IP.
 (May be suitable for use in other assays to be determined by the end user.)
Predicted Molecular Mass: 19.2kDa
Predicted isoelectric point: 5.6

Mouse CA12



Unique product Superb quality Client favorite Nicest service  ISO9001:2008;  ISO13485:2003; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnc.us | Http://www.uscnc.us
 Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnc.com | Http://www.uscnc.com

[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHSGSEF- ASDKSEGLAV LAVLIEIGSA NPSYDKIFSH LQHVKYKGQQ
VLIPGFNIEE LLPESPGEYY RYEGSLTTPP CYPTVLWTVF RNPVQISQEQ LLALETALYF
THMDDPTPRE MINNFRQVQK FDERLVYISF RQGLLTD TGL SLGIILSVAL AGVLGSI

Unique product Superb quality Client favorite Nicest service  ISO9001:2008;  ISO13485:2003; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us
Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com