



P93194Ra01 Collagen Type II Alpha 1 (COL2a1) Organism: Rattus norvegicus (Rat) *Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[DESCRIPTION]

Rat COL2a1 Protein Names: Collagen Type II Alpha 1 kDa Synonyms: COL2a1 94 66.2 Species: Rat 45 Size: 100ua 33 Source: Escherichia coli-derived 26 **Subcellular Location:** Secreted, extracellular space, extracellular matrix. 20 [PROPERTIES] Residues: Lys1226~Tyr1331 (Accession # P05539), 14.4 with N-terminal His-Tag. Grade & Purity: >95%, 15kDa as determined by SDS-PAGE reducing conditions. 15% SDS-PAGE Formulation: Supplied as lyophilized form in PBS, pH 7.4, containing 5% sucrose, 0.01% sarcosyl. **Endotoxin Level:** <1.0 EU per 1µg (determined by the LAL method). Applications: SDS-PAGE; WB; ELISA; IP. (May be suitable for use in other assays to be determined by the end user.) Predicted Molecular Mass: 13.5kDa Predicted isoelectric point: 7.1

Unique product Superb quality Client favorite Nicest service Q ISO9001:2008; Q ISO13485:2003; $\zeta \in$

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com





[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below. MGHHHHHHSGSEF- KSGDY WIDPNQGCTL DAMKVFCNME TGESCVYPNP ATVPRKNWWS SKSKEKKHIW FGETMNGGFH FSYGDGNLAP NTANVQMTFL RLLSTEGSQN ITYHCKNSIA Y

Unique product Superb quality Client favorite Nicest service @ ISO9001:2008; @ ISO13485:2003; \in