



# P91278Hu01 Fc Fragment Of IgG Low Affinity Illa Receptor (FcgR3A) Organism: Homo sapiens (Human) Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

5th Edition (Revised in January, 2013)

[ <u>DESCRIPTION</u> ] Protein Names: Fc Fragment Of IgG Low Affinity IIIa	Human FcgR3A	kDa
Receptor		94
Synonyms: FcgR3A, CD16A, FCG3, FCGR3, IGFR3		66.2
Species: Human		
Size: 100µg	-	45
Source: Escherichia coli-derived		
Subcellular Location: Cell membrane; Single-pass		33
type I membrane protein. Secreted.		26
[ PROPERTIES ]		20
Residues: Ser26~Ser108 (Accession # P08637), with		
N-terminal His-Tag.	-	20
Grade & Purity: >95%, 13kDa as determined by		
SDS-PAGE reducing conditions.		14.4
Formulation: Supplied as lyophilized form in PBS, pH		
7.4, containing 5% sucrose.		
<b>Endotoxin Level:</b> <1.0 EU per $1\mu g$ (determined by	15% SDS-PAG	ΞE
the LAL method).		
Applications: SDS-PAGE; WB; ELISA; IP.		
(May be suitable for use in other assays to be determined by the end user.)		
Predicted Molecular Mass: 10.8kDa		
Predicted isoelectric point: 5.9		

Unique product Superb quality Client favorite Nicest service (@ ISO9001:2008; @ ISO13485:2003;  $\zeta \in$ 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us Export Processing Zone, Wuhan 430056, P.R. China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com





# [PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

### [STORAGE AND STABILITY]

#### Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## [<u>SEQUENCES</u>]

The target protein is fused with N-terminal His-Tag, its sequence is listed below. MGHHHHHHSGSEF-S LVTCPLQCGI MWQLLLPTAL LLLVSAGMRT EDLPKAVVFL EPQWYRVLEK DSVTLKCQGA YSPEDNSTQW FHNESLISSQ AS

Unique product Superb quality Client favorite Nicest service @ ISO9001:2008; @ ISO13485:2003;  $\in$