



P90307Ra01 Galectin 7 (GAL7) Organism: Rattus norvegicus (Rat) *Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[DESCRIPTION]

Rat GAL7 Protein Names: Galectin 7 kDa 94 Synonyms: GAL7, Lgals7 66.2 Species: Rat 45 Size: 10ua 33 Source: Escherichia coli-derived 26 Subcellular Location: Cytoplasm, Nucleus, Secreted. 20 [PROPERTIES] Residues: Met1~Val128 (Accession # P97590), with 14.4 N-terminal His-Tag. Grade & Purity: >95%, 18kDa as determined by 15% SDS-PAGE SDS-PAGE reducing conditions. Formulation: Supplied as lyophilized form in PBS, pH 7.4, containing 5% sucrose, 0.01% sarcosyl. **Endotoxin Level:** <1.0 EU per 1µg (determined by the LAL method). Applications: SDS-PAGE; WB; ELISA; IP. (May be suitable for use in other assays to be determined by the end user.) Predicted Molecular Mass: 15.6kDa Predicted isoelectric point: 6.5

Unique product Superb quality Client favorite Nicest service Q ISO9001:2008; Q ISO13485:2003; $\zeta \in$

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com





[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[<u>SEQUENCES</u>]

The target protein is fused with N-terminal His-Tag, its sequence is listed below. MGHHHHHHSGS- MSATHHKTPL PQGVRLGTVM RIRGVVPDQA GRFHVNLLCG EEQEADAALH FNPRLDTSEV VFNTKQQGKW GREERGTGIP FQRGQPFEVL IITTEEGFKT VIGDDEYLHF HHRMPSSNVR SVEVGGDV

[<u>REFERENCES</u>]

- 1. Rondanino C., et al. (2011) Am. J. Physiol. Renal Physiol. 301:F622-33.
- 2. Saussez S., et al. (2006) Ann. Surg. Oncol. 13 (7): 999–1009.
- 3. Gerhard DS., et al. (2004) Genome Res. 14 (10B): 2121-7.
- 4. Kopitz J., et al. (2003) Oncogene 22 (40): 6277-88.

Unique product Superb quality Client favorite Nicest service 🙆 ISO9001:2008; 🧟 ISO13485:2003; $\zeta \in$