



P90050Mu01
Insulin Like Growth Factor 1 (IGF1)
Organism: Mus musculus (Mouse)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

1th Edition (Revised in February, 2012)

[DESCRIPTION]

Protein Names: Insulin Like Growth Factor 1

Gene Names: IGF1

Size: 100 μ g

Source: Recombinant

Expression Host: *E. coli*

Subcellular Location: Secreted

[PROPERTIES]

Residues: Gly33~Ala102 (Accession # Q8CAR0), with a N-terminal His-tag.

Grade & Purity: >97%, 9.2 kDa as determined by SDS-PAGE reducing conditions.

Form & Buffer: Supplied as lyophilized form in PBS, pH 7.4.

Endotoxin Level: <1.0 EU per 1 μ g(determined by the LAL method).

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 9.2 kDa

[PREPARATION]

Reconstitute in PBS.



[STORAGE AND STABILITY]

Storage: Store at 4°C for short term storage (1-2 weeks). Aliquot and store at -20°C or -80°C for long term storage. Avoid repeated freeze/thaw cycles.

Valid period: 12 months stored at -80°C.

[BACKGROUND]

The target protein is fused with a His-tag and its sequence is listed below. The first Met is an initiator amino acid. Moreover, Gly and Ser are added to improve the flexibility of N-terminus at both ends of the His-tag, which will increase the chelating ability of the tag to Ni-Sepharose during purification.

MGHHHHHSGSEF-GPETLCGA ELVDALQFVC GPRGFYFNKP TGYGSSIRRA PQTGIVDECC FRSCDLRRL
MYCAPLKPTK AA

[REFERENCES]

1. Rossetti,L., et al. (1996) J.Biol.Chem.271:203-208.
2. Carninci P., Hayashizaki Y.(1999) Methods Enzymol. 303:19-44.
3. The RIKEN Genome Exploration Research Group Phase II Team and the FANTOM Consortium (2001) Nature 409:685-690
4. Auernhammer,C.J.,et al. (1995) Euro.J.Endocrinol.133:635-645.
5. Delafontaine,P.,etal. (1996) J.Clin.Invest.97:139-145.
6. Carninci P., et al. (2000) Genome Res. 10:1617-1630
7. Rabkin,R., et al. (1995) J.Amer.Soc.Nephrol.6:1511-1518.

