

**P91073Mu02**  
**Interleukin 12 Receptor Beta 2 (IL12Rb2)**  
**Organism: Mus musculus (Mouse)**  
*Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

**[ DESCRIPTION ]**

**Protein Names:** Interleukin 12 Receptor Beta 2

**Synonyms:** IL12Rb2

**Species:** Mouse

**Size:** 100 $\mu$ g

**Source:** *Escherichia coli*-derived

**Subcellular Location:** Membrane; Single-pass type I membrane protein.

**[ PROPERTIES ]**

**Residues:** Thr455~Arg606 (Accession # P97378), with N-terminal His-Tag.

**Grade & Purity:** >95%, 21kDa as determined by SDS-PAGE reducing conditions.

**Formulation:** Supplied as lyophilized form in PBS, pH 7.4, containing 5% sucrose, 0.01% sarcosyl.

**Endotoxin Level:** <1.0 EU per 1 $\mu$ g (determined by the LAL method).


**Applications:** SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted Molecular Mass:** 19.1kDa

**Predicted isoelectric point:** 9.1

**Mouse IL12Rb2**



15% SDS-PAGE

*Unique product Superb quality Client favorite Nicest service*  ISO9001:2008;  ISO13485:2003; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us  
Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com

## **[ PREPARATION ]**

Reconstitute in sterile PBS, pH7.2-pH7.4.

## **[ STORAGE AND STABILITY ]**

**Storage: Avoid repeated freeze/thaw cycles.**

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## **[ SEQUENCES ]**

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHSGSEF- TWQPPK KADSAVREYI VEWRALQPGS ITKFPPHWLR IPPDNMSALI  
SENIKPYICY EIRVHALSES QGGCSSIRGD SKHKAPVSGP HITAITEKKE RLFISWTHPI  
FPEQRGCILH YRIYWKERDS TAQPELCEIQ YRRSQNSHPI SSLQPR