

P91134Ra02

Melanoma Associated Chondroitin Sulfate Proteoglycan (MCSP)

Organism: Rattus norvegicus (Rat)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

4th Edition (Revised in August, 2012)

[DESCRIPTION]

Protein Names: Melanoma Associated Chondroitin Sulfate Proteoglycan

Synonyms: MCSP, Ng2

Species: Rat

Size: 100µg

Source: *Escherichia coli*-derived

Subcellular Location: Cell membrane; Single-pass type I membrane protein; Extracellular side. Apical cell membrane; Extracellular side.

[PROPERTIES]

Residues: Gly1919~Lys2258 (Accession # Q00657), with two N-terminal Tags, His-tag and GST-tag.

Grade & Purity: >95%, 64 kDa as determined by SDS-PAGE reducing conditions.

Formulation: Supplied as liquid form in Phosphate buffered saline(PBS), pH 7.4.

Endotoxin Level: <1.0 EU per 1µg (determined by the LAL method).

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 63.6 kDa

Predicted isoelectric point: 6.6

[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYIADKHNLGGCPKERAIEISMLEGAVLDIRYGVSR IAYSKDFETLKVDFLSKLPEML
KMFEDRLCHKTYLNGDHVTHPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQG WQATFGGGDHPP
KSDGSTSGSGHHHHHSAGLVPRGSTAIGMKETA AAKFERQHMDSPDLGTLEVL FQGPLGSEF-GV FQLSMSDGAS PPIPMSLAVD
VLPSTIEVQL RAPLEVPQAL GRSSLSRQQL QVISDREEPD VAYRLTQG PL YGQVLVGGQP ASAFSQLQVD QGDVVFAFTN
FSSSQDHFV LALARGVNAS ATVNVTQAL LHVWAGGPWP QGTTLRDPT VLDASELANR TGSMRFRLL EGPRYGRVVR
VSQGRAESRT NQLVEDFTQQ DLEEGRLGLE VGRPEGRSTG PTGDRLTLEL QATGVPPAVA LLDFATEPYH AAKFYKV TLL
SVPEAARTET EKTGKSTPTG QPGQAASSPM PTVAKSGFLG FLEANMFSVI IPVCLVLLLL ALILPLLFYL RKRNK TGK

