

MAA447Po21

Monoclonal Antibody to C-Peptide (CP)

Organism Species: *Sus scrofa*; Porcine (Pig)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[PRODUCT INFORMATION]

Immunogen: C-Peptide-OVA

Clonality: Monoclonal

Clone number: E5

Host: Mouse

Immunoglobulin Type: IgG1 λ

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 500 μ g/mL

UOM: 200 μ g

[IMMUNOGEN INFORMATION]

Immunogen: Synthetic Peptide, C-Peptide conjugated to OVA.

Accession No.: CPA447Po21

[ANTIBODY SPECIFICITY]

The antibody is a mouse monoclonal antibody raised against C-Peptide. It has been selected for its ability to recognize C-Peptide in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.