

P90647Ra01
Mothers Against Decapentaplegic Homolog 1 (Smad1)
Organism: Rattus norvegicus (Rat)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[DESCRIPTION]

Protein Names: Mothers Against Decapentaplegic Homolog 1

Synonyms: Smad1, Mad1, Madh1

Species: Rat

Size: 50 μ g

Source: *Escherichia coli*-derived

Subcellular Location: Cytoplasm. Nucleus.

[PROPERTIES]

Residues: His273~Ser468 (Accession # P97588), with N-terminal His-Tag.

Grade & Purity: >95%, 25kDa as determined by SDS-PAGE reducing conditions.

Formulation: Supplied as lyophilized form in PBS, pH 7.4, containing 5% sucrose, 0.01% sarcosyl.

Endotoxin Level: <1.0 EU per 1 μ g (determined by the LAL method).


Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)




Predicted Molecular Mass: 23.8kDa

Predicted isoelectric point: 7.0

Rat Smad1



15% SDS-PAGE

Unique product Superb quality Client favorite Nicest service  ISO9001:2008;  ISO13485:2003; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnc.us | Http://www.uscnc.us
Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnc.com | Http://www.uscnc.com

[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHSGSEF-HWCSIVYY ELNNRVGERF HASSTSVLVD GFTDPSNNKN
RFCLGLLSNV NRNSTIENTR RHIGKGVHLY YVGGEVYAEC LSDSSIFVQS RNCNYHHGFH
PTTVCKIPSG CSLKIFNNQE FAQLLAQSVN HGFETEYELT KMCTIRMSFV KGWGAEYHRQ
DVTSTPCWIE IHLHGPLQWL DKVLTQMGSP HNPISVS

[REFERENCES]

1. Bucelli R.C., *et al.* (2008) J. Pharmacol. Exp. Ther. 324:1172-1180.
2. Pangas S.A., *et al.* (2008) Mol. Cell. Biol. 28:248-257.
3. Li P., *et al.* (2007) Nephrology (Carlton) 12:466-473.
4. Fuller M.L., *et al.* (2007) Ann. Neurol. 62:288-300.