CPA393Mi21 100µg OVA Conjugated Substance P (SP) Organism Species: Homo sapiens (Human) Mus musculus (Mouse) Rattus norvegicus (Rat) *Instruction manual* 

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

# Coud-Clone Corp.

#### [PROPERTIES]

Source: Protein Conjugation Original Peptide Residues: Arg58-Met68 Purity: >90% Traits: Freeze-dried powder. Buffer formulation: PBS, pH7.4. Applications: Immunogen; Coating Antigen; Blocking Peptide; ELISA.

### [ <u>USAGE</u> ]

Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

#### [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.


**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.


## [<u>SEQUENCE</u>]

RPK PQQFFGLM



#### [IDENTIFICATION]





1304 Langham Creek Dr. Suite 226, Huustun, TX 77084, USA | 001-888-960-7402 | www.cloud-cloue.us | mail@cloud-cloue.as Export Processing Zone, Wuhan, Hubei 430056, PRC | 0086-008-880-0687 | www.cloud-clone.com | mail@cloud-clone.com