



P95757Hu02 Peroxiredoxin 2 (PRDX2) Organism: Homo sapiens (Human) *Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[DESCRIPTION]

Protein Names: Peroxiredoxin 2 Human PRDX2 Synonyms: PRDX2, NKEFB, TDPX1 kDa 94 66 2 Species: Human 45 Size: 100ua Source: Escherichia coli-derived 33 Subcellular Location: Cytoplasm. 26 [PROPERTIES] 20 **Residues:** Ala6~Tyr164 (Accession # P32119), with two N-terminal Tags, His-tag and GST-tag. 14.4 Grade & Purity: >95%, 46kDa as determined by SDS-PAGE reducing conditions. **Formulation:** Supplied as lyophilized form in PBS, pH 15% SDS-PAGE 7.4, containing 5% sucrose, 0.01% sarcosyl. **Endotoxin Level:** <1.0 EU per 1µg (determined by the LAL method). Applications: SDS-PAGE; WB; ELISA; IP. (May be suitable for use in other assays to be determined by the end user.) Predicted Molecular Mass: 44.7kDa Predicted isoelectric point: 6.5







[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[<u>SEQUENCES</u>]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MRNKKFELGL EFPNLPYYID GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV DFLSKLPEML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK KRIEAIPQID KYLKSSKYIA WPLQGWQATF GGGDHPPKSD GSTSGSGHHH HHHSAGLVPR GSTAIGMKET AAAKFERQHM DSPDLGTLEV LFQ GPLGSEF- ARIGK PAPDFKATAV VDGAFKEVKL SDYKGKYVVL FFYPLDFTFV CPTEIIAFSN RAEDFRKLGC EVLGVSVDSQ FTHLAWINTP RKEGGLGPLN IPLLADVTRR LSEDYGVLKT DEGIAYRGLF IIDGKGVLRQ ITVNDLPVGR SVDEALRLVQ AFQY

Unique product Superb quality Client favorite Nicest service 🙆 ISO9001:2008; 🧟 ISO13485:2003; $\zeta \in$