

PAA579Hu01

Polyclonal Antibody to Aquaporin 1, Colton Blood Group (AQP1)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PRODUCT INFORMATION]

Immunogen: AQP1, Human **Purification:** Affinity Chromatography.

Clonality: Polyclonal Applications: WB, ICC, IHC-P, IHC-F, ELISA

Host: Rabbit Concentration: 200µg/mL

Immunoglobulin Type: IgG **UOM**: 100μg

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant AQP1 (Val201~Lys269) with two N-terminal Tags,

His-tag and GST-tag expressed in E.coli.

Accession No.: RPA579Hu01

[ANTIBODY SPECIFITY]

The antibody is a rabbit polyclonal antibody raised against AQP1. It has been selected for its ability to recognize AQP1 in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:50-400

Immunocytochemistry in formalin fixed cells: 1:50-500

Immunohistochemistry in formalin fixed frozen section: 1:50-500

Immunohistochemistry in paraffin section: 1:10-100 Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.



[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[QUALITY CONTROL]

Content: The quality control contains recombinant AQP1 (Val201~Lys269) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemilumescent (ECL).


Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading Buffer: 100mM Tris(pH8.8), 2% SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, NaN₃ 0.02%.


[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.

[IMAGES]



Used in Western Blot, Sample:
Recombinant AQP1, Human




Used in DAB staining on fromalin fixed paraffin- embedded kidney tissue




PAA579Hu01

Polyclonal Antibody to Aquaporin 1, Colton Blood Group (AQP1) Organism Species: Homo sapiens (Human) Instruction manual


FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES




Used in DAB staining on fromalin fixed paraffin- embedded liver tissue



Used in DAB staining on fromalin fixed paraffin- embedded liver cancer tissue



Used in DAB staining on fromalin fixed paraffinembedded prostate glands cancer tissue



Used in DAB staining on fromalin fixed paraffin- embedded glioma tissue




PAA579Hu01

Polyclonal Antibody to Aquaporin 1, Colton Blood Group (AQP1)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES



Used in DAB staining on fromalin fixed paraffinembedded rectum cancer tissue