

PAA284Hu01

Polyclonal Antibody to Heparin Cofactor II (HCII)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PRODUCT INFORMATION]

Immunogen: HCII, Human Purification: Affinity Chromatography.

Clonality: Polyclonal Applications: WB, ICC, IHC-P, IHC-F, ELISA

Host: Rabbit Concentration: 200µg/mL

Immunoglobulin Type: IgG **UOM**: 100μg

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant HCII (Asp276~Ser499) with two N-terminal Tags,

His-tag and T7-tag expressed in E.coli.

Accession No.: RPA284Hu01

[ANTIBODY SPECIFITY]

The antibody is a rabbit polyclonal antibody raised against HCII. It has been selected for its ability to recognize HCII in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:50-400

Immunocytochemistry in formalin fixed cells: 1:50-500

Immunohistochemistry in formalin fixed frozen section: 1:50-500

Immunohistochemistry in paraffin section: 1:10-100 Enzyme-linked Immunosorbent Assay: 1:100-5000

Optimal working dilutions must be determined by end user.




[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.


[STORAGE]

Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.


[IMAGES]



Used in Western Blot, Sample: Recombinant HCII, Human



Used in DAB staining on fromalin fixed paraffin- embedded liver tissue



Used in Western Blot, Sample:

Lane1: Human Serum

Lane2: Porcine Lung Tissue

Lane3: Porcine Liver Tissue

Lane4: Porcine Kidney Tissue

Primary Ab: 1µg/mL Rabbit Anti-Human HCII Ab

Second Ab: 1:1000 Dilution of HRP-Linked Guinea pig

Anti-Rabbit Ab (Catalog: SECab05)