

RPL231Mu01 10µg

Recombinant Alpha-2-Glycoprotein 1, Zinc Binding (aZGP1)

**Organism Species: Mus musculus (Mouse)** 

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

11th Edition (Revised in May, 2016)

#### [PROPERTIES]

**Source:** Prokaryotic expression.

Host: E. coli

Residues: Glu178~Gln307 Tags: N-terminal His-Tag

Tissue Specificity: Liver, Kidney, Plasm.

Subcellular Location: Secreted.

**Purity: >95%** 

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 1mM DTT, 5% trehalose, 0.01%

sarcosyl and Proclin300.

Original Concentration: 200ug/mL

Applications: SDS-PAGE; WB; ELISA; IP; CoIP; Reporter Assays; Purification;

Amine Reactive Labeling.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 7.1

Predicted Molecular Mass: 16.6kDa

Accurate Molecular Mass: 17kDa as determined by SDS-PAGE reducing conditions.

#### [USAGE]

Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.



### [STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

# [SEQUENCE]

PPTVTITSRV IPGGNRIFKC LAYGFYPQRI SLHWNKANKK LAFEPERGVF PNGNGTYLSW AEVEVSPQDI DPFFCLIDHR GFSQSLSVQW DRTRKVKDEN NVVAQPQ

## [ IDENTIFICATION ]




Figure 1. SDS-PAGE