

RPG442Hu01 100µg

Recombinant Cysteine Conjugate Beta Lyase 2 (CCbL2)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Met1~Ser454

Tags: N-terminal His Tag

Subcellular Location: Mitochondrion, Cytoplasm

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.4

Predicted Molecular Mass: 55.1kDa

Accurate Molecular Mass: 52kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MFLAQRSLCS LSGRAKFLKT ISSSKILGFS TSAKMSLKFT NAKRIEGLDS NWWIEFTKL AADPSVVNLG QGFPDISPPT YVKEELSKIA
AIDSLNQYTR GFGHPSLVK ALSYLYEKLY QKQIDSNKEI LVTVGAYGSL FNTIQALIDE GDEVILIVP FYDCYEPMVR MAGATPVFIP
LRSKPVYGKR WSSSDWTLDP QELESKFNS KTKAILLNTF HNPLGKVYNR BELQVIADLC IKYDTLCISD EYEWLVYS GNKHLKIATF
PGMWERTITI GSAGKTFSVT GWKLGWSIGP NHLIKHLQT VQQNTIYTCA TPLQEALQA FWIDIKRMDD PECYFNSLPK ELEVKRDRM
VRLLESVGLK PIVPDGGYFI IADVSLDPD LSDMKNNPEY DYKFKWMT KHKLSAIPV SAFCNSETKS QFEKVRFCF IKKDSTLDA
EEIKAWSV QKS

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.