

RPA303Rb01 100µg

Recombinant Galectin 3 (GAL3)

Organism Species: Oryctolagus cuniculus (Rabbit)

Instruction manual

kDa

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Ala2~lle242

Tags: Two N-terminal Tags, His-tag and GST-tag

Accession: P47845

Host: E. coli

Subcellular Location: Cytoplasm. Nucleus.

Secreted.

Purity: >90%

Endotoxin Level: <1.0EU per 1μg (determined by the LAL method).

Formulation: Supplied as lyophilized form in PBS,

pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 9.0

Predicted Molecular Mass: 55.4kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

70 44 33 26 22 18 14 10

[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

ADGFSLNDA LSGSGHPPNQ GWPGPWGNQP AGPGGYPGAA YPGAYPGHAP GAYPGQAPPG PYPGPGAHGA YPGQPGGPGA YPSPGQPSGA GAYPGASPYS ASAGPLPVPY DLPLPGGVMP RMLITIVGTV KPNANRLALD FKRGNDVAFH FNPRFNENNR RVIVCNTKVD NNWGREERQT TFPFEIGKPF KIQVLVEPDH FKVAVNDAHL LQYNHRMRNL KEINKLGISG DIQLTSASHA MI