

RPG481Mu01 100ug

Recombinant Growth Arrest Specific Protein 8 (GAS8)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Met1~Thr478

Tags: N-terminal His Tag

Subcellular Location: Secreted

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 7.7

Predicted Molecular Mass: 60kDa

Accurate Molecular Mass: 60kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.


Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MAPKKKGKKG KAKGTAIVDG VAPEDMTKEQ VEEHVARIRE ELDREREERN
YFQLERDKI HTFWEITRRQ LEEKKAELRN KDREMEEAEE RHQVEIKVYK Q
KVKHLLYE HQNNLAEVKA EGTVVMKLAQ KEHRTQEGAL RKDMRVLKVE LK
EQELANE VVIKNLCLKQ AEEITKMRND FERQVREIEA KYDKKMKMLR DEL
DLRRKT EIHEVEERKN GQISTLMQRH EEAFTDIKNY YNDITLNNLA LINS
LKEQM EDMRKKEEHM EREMAEVTLQ NRRLADPLQK AKDEMNMOMQK RLGNH
ERDK QILVCTKARL KVAERELKDL KWEHEVLEQR FIKVQQEREE LYRKFA
DAI QEVQQKTGFK NLLLERKLQA LNAAVEKREV QFNEVLAASN LDPTALT
LV SRKLEDVLES KNTTIKDLQY ELARVCKAHN DLLRTYEAKL LAFGIPLD
N VGFKPLETAV IGQTLGQGPA GLVGAPT

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.