

**RPC252Hu01 10µg**

**Recombinant Heterogeneous Nuclear Ribonucleoprotein A1 (HNRNPA1)**

**Organism Species: *Homo sapiens (Human)***

***Instruction manual***

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

---

13th Edition (Revised in Aug, 2023)

**[ PROPERTIES ]**

**Source:** Prokaryotic expression

**Host:** *E.coli*

**Residues:** Ser2~Phe372

**Tags:** N-terminal His Tag

**Subcellular Location:** Nucleus, Cytoplasm

**Purity:** > 97%

**Traits:** Freeze-dried powder

**Buffer formulation:** 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% skl, 5%Trehalose.

**Original Concentration:** 50µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 9.2

**Predicted Molecular Mass:** 42.3kDa

**Accurate Molecular Mass:** 43kDa as determined by SDS-PAGE reducing conditions.

**[ USAGE ]**

Reconstitute in ddH<sub>2</sub>O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

**[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

**[ SEQUENCE ]**

SKSESPKEP EQLRKLFIGG LSFETTDESL RSHFEQWGL TDCVVMRDPN  
TKRSRGFGFV TYATVEEVDA AMNARPHKVD GRVVEPKRAV SREDSQRPGA  
HLTVKKIFVG GIKEDTEHH LRDYFEQYGK IEVIEIMTDR GSGKKRGFAF  
VTFDDHDSVD KIVIQKYHTV NGHNCVVRKA LSKQEMASAS SSQRGRSGSG  
NFGGGRGGGF GGNDNFGRGG NFSGRGGFGG SRGGGGYGGG GDGYNGFGND  
GGYGGGGPGY SGGSRGYGSG GQGYGNQGS YGGSGSYDSY NNGGGGGFGG  
GSGSNFGGGG SYNDFGNYYN QSSNFGPMKG GNFGGRSSGP YGGGGQYFAK  
PRNQGGYGGG SSSSSYGSGR RF

**[ IDENTIFICATION ]**




Figure. SDS-PAGE

**[ IMPORTANT NOTE ]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.