

RPE768Hu01 100µg

Recombinant Interleukin 20 Receptor Alpha (IL20Ra)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES.

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Met269~Trp545

Tags: Two N-terminal Tags, His-tag and GST-tag

Accession: Q9UHF4

Host: E. coli

Subcellular Location: Membrane; Single-pass

type I membrane protein.

Purity: >90%

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).

Formulation: Supplied as lyophilized form in PBS,


pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 4.1

Predicted Molecular Mass: 61.1kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

MG YSIYRYIHVG KEKHPANLIL IYGNEFDKRF FVPAEKIVIN FITLNISDDS KISHQDMSLL GKSSDVSSLN DPQPSGNLRP PQEEEEVKHL GYASHLMEIF CDSEENTEGT SLTQQESLSR TIPPDKTVIE YEYDVRTTDI CAGPEEQELS LQEEVSTQGT LLESQAALAV LGPQTLQYSY TPQLQDLDPL AQEHTDSEEG PEEEPSTTLV DWDPQTGRLC IPSLSSFDQDSEGCEPSEGD GLGEEGLLSR LYEEPAPDRP PGENETYLMQ FMEEW