

RPK166Mu01 100µg

Recombinant Leucine Rich Repeat LGI Family, Member 3 (LGI3)

Organism Species: Mus musculus (Mouse)

Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)



# [PROPERTIES]

Source: Prokaryotic expression

Host: E.coli

Residues: Lys261~Ala548

Tags: N-terminal His Tag

Subcellular Location: Secreted

**Purity:** > 97%

Traits: Freeze-dried powder

**Buffer formulation:** 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.9

Predicted Molecular Mass: 36.8kDa

**Accurate Molecular Mass:** 34kDa as determined by SDS-PAGE reducing conditions.

#### [USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

#### [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

# [ SEQUENCE ]



```
KWDYVERQLR DYDRIPAPSA VHCKPMVVDG QLYVVVAQLF
GGSYIYHWDP NTTRFTKLQD IDPQRVRKPN DLEAFRIDGD WFFAVADSSK
AGATSLYRWH QNGFYSHQAL HAWHRDTDLE FVDGEGKPRL IVSSSSQAPV
IYQWSRSQKQ FVAQGEVTQV PDAQAVKHFR AGRDSYLCLS RYIGDSKILR
WEGTRFSEVQ ALPSRGSLAL QPFLVGGHRY LALGSDFSFT QIYQWDEGRQ
KFVRFQELAV QAPRAFCYMP AGDAQLLLAP SFKGQTLVYR HVVVDLSA
```

# [ IDENTIFICATION ]




Figure. SDS-PAGE

### [ IMPORTANT NOTE ]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.