

RPB853Hu01 50µg Recombinant Matrix Metalloproteinase 19 (MMP19) Organism Species: *Homo sapiens (Human) Instruction manual* 

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

# Cond-Clone Corp.

### [PROPERTIES]

Source: Prokaryotic expression

Host: E.coli

**Residues:** Tyr98~Tyr508

Tags: N-terminal His Tag

Subcellular Location: Secreted, Extracellular matrix

**Purity:** > 95%

Traits: Freeze-dried powder

Buffer formulation: 100mMNaHCO<sub>3</sub>, 500mMNaCl, pH8.3, containing 0.01% SKL, 5%

Trehalose.

Original Concentration: 300µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.3

Predicted Molecular Mass: 50.2kDa

Accurate Molecular Mass: 50kDa as determined by SDS-PAGE reducing conditions.

#### [<u>USAGE</u>]

Reconstitute in  $ddH_2O$  to a concentration of 0.1-1.0 mg/mL. Do not vortex.

#### [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.


#### [<u>SEQUENCE</u>]

## Cloud-Clone Corp.

YLL

LGRWRKKHLT FRILNLPSTL PPHTARAALR QAFQDWSNVA PLTFQEVQAG AADIRLSFHG RQSSYCSNTF DGPGRVLAHA DIPELGSVHF DEDEFWTEGT YRGVNLRIIA AHEVGHALGL GHSRYSQALM APVYEGYRPH FKLHPDDVAG IQALYGKKSP VIRDEEEET ELPTVPPVPT EPSPMPDPCS SELDAMMLGP RGKTYAFKGD YVWTVSDSGP GPLFRVSALW EGLPGNLDAA VYSPRTQWIH FFKGDKVWRY INFKMSPGFP KKLNRVEPNL DAALYWPLNQ KVFLFKGSGY WQWDELARTD FSSYPKPIKG LFTGVPNQPS AAMSWQDGRV YFFKGKVYWR LNQQLRVEKG YPRNISHNWM HCRPRTIDTT PSGGNTTPSG TGITLDTTLS ATETTFEY

#### [IDENTIFICATION]



#### [IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.