

RPA754Hu01 50µg

#### Recombinant Myosin Heavy Chain 1, Skeletal Muscle, Adult (MYH1) Organism Species: Homo sapiens (Human)

#### Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

11th Edition (Revised in May, 2016)

# [PROPERTIES]

Source: Prokaryotic expression. Host: E. coli Residues: Tyr585~Ala794 Tags: N-terminal His-Tag Tissue Specificity: Skeletal Muscle. Subcellular Location: Cytoplasm. **Purity: >98%** Traits: Freeze-dried powder Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% sarcosyl, 5% Trehalose and Proclin300. Original Concentration: 200ug/mL Applications: SDS-PAGE; WB; ELISA; IP; CoIP; Purification; Amine Reactive Labeling. (May be suitable for use in other assays to be determined by the end user.) Predicted isoelectric point: 9.6 Predicted Molecular Mass: 27.0kDa Accurate Molecular Mass: 27kDa as determined by SDS-PAGE reducing conditions.

### [ <u>USAGE ]</u>

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.



### [STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

## [ <u>SEQUENCE</u> ]

YAGTVD YNIAGWLDKN KDPLNETVVG LYQKSAMKTL ALLFVGATGA EAEAGGGKKG GKKKGSSFQT VSALFRENLN KLMTNLRSTH PHFVRCIIPN ETKTPGAMEH ELVLHQLRCN GVLEGIRICR KGFPSRILYA DFKQRYKVLN ASAIPEGQFI DSKKASEKLL GSIDIDHTQY KFGHTKVFFK AGLLGLLEEM RDEKLAQLIT RTQA

#### [IDENTIFICATION]

and the second second	
	kDa 70
	44
100	33
-	26
	22
	18
	14
	10

Figure 1. SDS-PAGE