

RPC728Hu01 50µg

Recombinant Phospholipid Transfer Protein (PLTP)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Val164~Val493

Tags: N-terminal His Tag

Subcellular Location: Nucleus, Secreted

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% skl, 5%Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.5

Predicted Molecular Mass: 40.2kDa

Accurate Molecular Mass: 40kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

VYDFLST FITSGMRFL NQICPVLYH AGTVLLNSLL
 DTPVVRSSVD ELVGIDYSLM KDPVASTSNL DMDFRGAFFP LTERNWSLPN
 RAVEPQLQEE ERMVYVAFSE FFFDSAMESY FRAGALQLLL VGDKVPDLD
 MLLRATYFGS IVLLSPAVID SPLKLELRVL APPRCTIKPS GTTISVTASV
 TIALVPPDQP EVQLSSMTMD ARLSAKMALR GKALRTQLDL RRFRIYSNHS
 ALESLALIPL QAPLKTMLQI GVMPMLNERT WRGVQIPLPE GINFEVHEVVT
 NHAGFLTIGA DLHFAKGLRE VIEKNRPADV RASTAPTST AAV

[IDENTIFICATION]




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.