

RPD935Mu01 100ug

Recombinant Proteasome Assembly Chaperone 2 (PSMG2)

Organism Species: Mus musculus (Mouse)

Instruction manual

kDa

70

44

33

26

22

18

14

10

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES.

10th Edition (Revised in Jan. 2014)

[PROPERTIES]

Residues: Met1~Phe264

Tags: Two N-terminal Tags, His-tag and T7-tag

Accession: Q9EST4

Host: E. coli

Subcellular Location: Nucleus.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).

Formulation: Supplied as lyophilized form in 20mM Tris,

150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT,

Predicted Molecular Mass: 33.2kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

0.01% sarcosyl, 5% trehalose, and preservative. 15% SDS-PAGE Predicted isoelectric point: 5.8

[USAGE]

Reconstitute in sterile ddH₂O.



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

MFVPCGESVP DLTNFTLLMP AVSVGNVGQL AIDLIISTLN MCKIGYFYTD CLVPMVGNNP YATEEENSNE LSINTEVYSL PSKKLVVLQL RSIFIKYKSK SFCEKLLAWV ESSGCARIIV LSSSHSYHRN DAQLRSTPFR YLLTPCLQKS VQNKIKSLNW LEMEKSRCIP EMSDSEFCIR IPGGGITKTL YDESCSKEIQ MAVLLKFVSE GDNIPDAVSL VEYLNEWLQI IKPCNDGPMA SALPWKIPSS WRLLFGSGLP PALF