

RPB786Ra01 10μg Recombinant Receptor Interacting Serine Threonine Kinase 2 (RIPK2) Organism Species: *Rattus norvegicus (Rat) Instruction manual*

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

Cloud-Clone Corp.

[PROPERTIES]

Source: Prokaryotic expression Host: *E.coli* Residues: Met1~Leu294 Tags: N-terminal His Tag Subcellular Location: Cytoplasm Purity: > 97% Traits: Freeze-dried powder Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose . Original Concentration: 50µg/mL Applications: Positive Control; Immunogen; SDS-PAGE; WB. (May be suitable for use in other assays to be determined by the end user.) Predicted isoelectric point: 6.5 Predicted Molecular Mass: 37.5kDa

Accurate Molecular Mass: 39kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH_2O to a concentration of 0.1-0.25 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.


Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]



MNGDAICSAL PHIPYHKLAD LHYLSRGASG TVSSARHADW RVRVAVKHLH IHTPLLDSER NDILREAEIL HKARFSYILP ILGICNEPEF LGIVTEYMPN GSLNELLHRK TEYPEIAWPL RFRILHEIAL GVNYLHNMNP PLLHHDLKTQ NILLDNEFHV KIADFGLSKW RMMSLSQSRS YKSAPEGGTI IYMPPENYEP GQKSRASVKH DIYSYAVIMW EVLSRKQPFE EVTNPLQIMY SVSQGHRPNT SEENLPFDIP HRGLMISLIQ SGWAQNPDER PSFLKCLIEL EPVL

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.