

RPA116Po01 100µg

Recombinant Regulated On Activation

In Normal T-Cell Expressed And Secreted (RANTES)

Organism Species: scrofa; Porcine (Pig)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Ser24~Ser91

Tags: Two N-terminal Tags, His-tag and T7-tag

Accession: Q2EN89

Host: E. coli

Subcellular Location: Cytoplasm.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).

Formulation: Supplied as lyophilized form in PBS,

pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 9.1

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 11.6kDa Applications: SDS-PAGE; WB; ELISA; IP.

[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

SPYASDT TPCCFSYLSR PLPRAHLQEY FYTSSKCSMA AVVFITRKNR QVCANPEKKW VREYINSLEL S