

RPG850Ra01 100μg Recombinant Secretagogin (SCGN) Organism Species: *Rattus norvegicus (Rat) Instruction manual*

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

Coud-Clone Corp.

[PROPERTIES]

Source: Prokaryotic expression Host: *E.coli* Residues: Met1~Pro276 (Accession # Q6R556) Tags: N-terminal His Tag Subcellular Location: Secreted, Cytoplasm Purity: > 97% Traits: Freeze-dried powder Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose. Original Concentration: 200µg/mL Applications: Positive Control; Immunogen; SDS-PAGE; WB. (May be suitable for use in other assays to be determined by the end user.) Predicted isoelectric point: 5.2 Predicted Molecular Mass: 36.0kDa Accurate Molecular Mass: 34kDa as determined by SDS-PAGE reducing conditions.

[<u>USAGE</u>]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.


Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]



MDNAHRQTQA	HLDAACFWQI	WQRFDKDEKG	YIKETELDAF	FDDLLAKFGI
EDTLMEENVQ	KMKEQLMVGH	DISKEGRILM	KELASMFLSE	DENFLLFFRL
ETPLDNSVEF	MQIWRKYDAD	SSGFISAAEL	SNFLRDLFLH	HKKVI SEAEL
EEYTSTMEKI	FDRNKDGRLD	LNDLARILAL	QENFLLQFKM	DASSTEERKR
DFEKIFAHYD	VSKTGALEGP	EVDGFVKDMM	ELVQPSISGV	DLDKFREILL
RHCDVNKDGK	IQKSELALCL	GLKINP		

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.