#### RPC796Hu01 100µg Recombinant Tumor Necrosis Factor Alpha Induced Protein 3 (TNFaIP3) Organism Species: Homo sapiens (Human)

#### Instruction manual

#### FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

### [PROPERTIES]


Residues: Met1~Glu300 Tags: Two N-terminal Tags, His-tag and T7-tag Accession: P21580 Host: *E. coli* Subcellular Location: Cytoplasm, Lysosome, Nucleus. Purity: >95% Endotoxin Level: <1.0EU per 1µg (determined by the LAL method). Formulation: Supplied as Iyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl. Predicted isoelectric point: 8.6 Predicted Molecular Mass: 38.9kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

# [<u>USAGE</u>]

Reconstitute in sterile PBS, pH7.2-pH7.4.





### [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

# [SEQUENCES]

The sequence of the target protein is listed below.

MAEQVLPQAL YLSNMRKAVK IRERTPEDIF KPTNGIIHHF KTMHRYTLEM FRTCQFCPQF REIIHKALID RNIQATLESQ KKLNWCREVR KLVALKTNGD GNCLMHATSQ YMWGVQDTDL VLRKALFSTL KETDTRNFKF RWQLESLKSQ EFVETGLCYD TRNWNDEWDN LIKMASTDTP MARSGLQYNS LEEIHIFVLC NILRRPIIVI SDKMLRSLES GSNFAPLKVG GIYLPLHWPA QECYRYPIVL GYDSHHFVPL VTLKDSGPEI RAVPLVNRDR GRFEDLKVHF LTDPENEMKE