

RPB040Hu01 50µg

Recombinant Vimentin (VIM)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Ser2~Glu466

Tags: N-terminal His Tag

Subcellular Location: Cytoplasm

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% Sarcosyl, 5% Trehalose.

Original Concentration: 500µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.2

Predicted Molecular Mass: 55.0kDa

Accurate Molecular Mass: 52kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

STRSVSSSS	YRRMFGGPGT	ASRPSSSRSY	VTTSTRTYSL	GSALRPSTSR
SLYASSPGGV	YATRSSAVRL	RSSVPGVRL	QDSVDFSLAD	AINTEFKNTR
TNEKVELQEL	NDRFANYIDK	VRFLEQQNKI	LLAELEQLKG	QGKSRLGDLY
EEEMRELRRQ	VDQLTNDKAR	VEVERDNLAE	DIMRLREKLQ	EEMLQREEAE
NTLQSFRQDV	DNASLARLDL	ERKVESLQEE	IAFLKKLHEE	EIQELQAQIQ
EQHVQIDVDV	SKPDLTAALR	DVRQQYESVA	AKNLQEAEW	YKSKFADLSE
AANRNNDALR	QAKQESTEYR	RQVQSLTCEV	DALKGTNESL	ERQMREMEEN
FAVEAANYQD	TIGRLQDEIQ	NMKEEMARHL	REYQDLLNVK	MALDIEIATY
RKLLEGEESR	ISLPLPNFSS	LNLRETNLDS	LPLVDTHSKR	TLLIKTVETR
DGQVINETSQ	HHDDLE			

[IDENTIFICATION]




Figure . Gene Sequencing (extract)



Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.