ු <u>Cloud-Clone Corp</u>.

RPA394Rb01 100µg Recombinant Tissue Factor Pathway Inhibitor (TFPI) Organism Species: Oryctolagus cuniculus (Rabbit) *Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

[PROPERTIES]

Residues: Lys39~Gly275 (Accession # P19761), with two N-terminal Tags, His-tag and GST-tag. Host: *E. coli* Subcellular Location: Secreted. Purity: >95% Endotoxin Level: <1.0EU per 1µg (determined by the LAL method). Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 1mM DTT, 5% trehalose, 0.05% sarcosyl and preservative. Predicted isoelectric point: 7.9 Predicted Molecular Mass: 59.2kDa

(May be suitable for use in other assays to be determined by the end user.)


[<u>USAGE</u>]

Reconstitute in sterile PBS, pH7.2-pH7.4.

Applications: SDS-PAGE; WB; ELISA; IP.

kDa

9th Edition (Revised in Jul, 2013)



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[<u>SEQUENCES</u>]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MSPILGYWKI KGLVQPTRLL LEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV DFLSKLPEML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK KRIEAIPQID KYLKSSKYIA WPLQGWQATF GGGDHPPKSD GSTSGSGHHH HHHSAGLVPR GSTAIGMKET AAAKFERQHM DSPDLGTLEV LFQGPLGSEF-KP PLQKPTHSFC AMKVDDGPCR AYIKRFFFNI LTHQCEEFIY GGCEGNENRF ESLEECKEKC ARDYPKMTTK LTFQKGKPDF CFLEEDPGIC RGYITRYFYN NQSKQCERFK YGGCLGNLNN FESLEECKNT CENPTSDFQV DDHRTQLNTV NNTLINQPTK APRRWAFHGP SWCLPPADRG LCQANEIRFF YNAIIGKCRP FKYSGCGGNE NNFTSKKACI TACKKGFIPK SIKGG

[REFERENCES]

- 1. Wesselschmidt R.L., et al. (1990) Nucleic Acids Res. 18:6440-6440.
- 2. Warn-Cramer B.J., et al. (1992) Nucleic Acids Res. 20:3548-3548.
- 3. Belaaouaj A., et al. (1993) Thromb. Res. 69:547-553.
- 4. Luo J.G., et al. (2009) Zhonghua Xin Xue Guan Bing Za Zhi 37:1113-1118.

Designed by Cloud-Clone Corp., Assembled by Uscn Life Science Inc. ISO9001:2008; ISO13485:2003 11271 Richmond Avenue Suite H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Http://www.cloud-clone.us | E-mail: mail@cloud-clone.us
Export Processing Zone Building F. Wuhan, Hubel 430056, PRC | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-8551 | Http://www.uscnk.com | E-mail: mail@cloud-clone.us