

PAA166Ra01**Polyclonal Antibody to Pepsinogen C (PGC)****Organism Species: Rattus norvegicus (Rat)*****Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[PRODUCT INFORMATION]**Immunogen:** PGC, Rat**Clonality:** Polyclonal**Host:** Rabbit**Immunoglobulin Type:** IgG**Purification:** Affinity Chromatography.**Applications:** WB, ICC, IHC-P, IHC-F, ELISA**Concentration:** 200µg/mL**UOM:** 100µg**[IMMUNOGEN INFORMATION]****Immunogen:** Recombinant PGC (Ser17~Val392) expressed in *E.coli*.**USCN Accession No.:** RPA166Ra01**Sequence:** The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHSGSEF-SLLR VPLRKMKSIR ETMKEQGVLK DFLKTHKYDP GQKYHFGNFG
DYSVLYEPMA YMDASYFGEI SIGTPPQNFL VLFDTGSSNL WVSSVYCQSE ACTTHARFNP
SKSSTYYTEG QTFSLQYGTG SLTGFFGYDT LTVQSIQVPN QEFGLSENEP GTNFVYAQFD
GIMGLAYPGL SSGGATTALQ GMLGEGALSQ PLFGVYLGSQ QGSNGGQIVF GGVDKNLYTG
EITWVPVTQE LYWQITIDDF LIGDQASGWC SSQGCQGIVD TGTSLLVMPA QYLSELLQTI
GAQEGEYGEY FVSCDSVSSL PTLSFVLNGV QFPLSPSSYI IQEDNFCMVG LESISLTSES
GQPLWILGDV FLRSYYAIFD MGNKVG/LAT SV

[ANTIBODY SPECIFICITY]

The antibody is a rabbit polyclonal antibody raised against PGC. It has been selected for its ability to recognize PGC in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500


Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200


Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[**IMAGES**]



Used in DAB staining on formalin fixed
paraffin-embedded kidney tissue



Used in Western Blot, Sample:
Recombinant PGC, Rat

[**CONTENTS**]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN_3 , 50% glycerol.

[**QUALITY CONTROL**]

Content: The quality control contains recombinant PGC (Ser17~Val392) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.
5uL per well when used in enhanced chemiluminescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading Buffer: 100mM Tris(pH8.8), 2% SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, NaN_3 0.02%.

[**STORAGE**]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.