Coud-Clone Corp.

SAA544Eq09 HRP-Linked Rabbit Anti-Horse IgG Antibody Organism Species: Equus caballus; Equine (Horse) *Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

11th Edition (Revised in May, 2016)

[PROPERTIES]

Source: Antibody labeling Host: Rabbit Purification: Antigen-specific Affinity Chromatography. Label: HRP Traits: Liquid Concentration: 200µg/mL UOM: 100µg Applications: WB; IP; IHC-P; IHC-F; ELISA.

[IMMUNOGEN]

Immunogen: IgG, Equine

[APPLICATIONS]

Western blotting: 1:2000-10000 Immunoprecipitation: 1:2000-10000 Immunohistochemistry in formalin fixed frozen section: 1:200-1000 Immunohistochemistry in paraffin section: 1:200-1000 Enzyme-linked Immunosorbent Assay: 1:4000-15000 Optimal working dilutions must be determined by end user.

[FORMULATION]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

Coud-Clone Corp.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 4°C for frequent use.

Aliquot and store at -20°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[IDENTIFICATION]




Figure 1: Used in Western Blot, Sample: Equine Serum

Ab: 1:2000 Dilution of HRP-Linked Rabbit Anti-Horse IgG Ab (Catalog: SAA544Eq09)