Coud-Clone Corp.

CPA441Hu21 100µg OVA Conjugated Luteinizing Hormone (LH) Organism Species: Homo sapiens (Human) *Instruction manual* 

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

11th Edition (Revised in May, 2016)

#### [PROPERTIES]

Source: Protein Conjugation Original Peptide Residues: Lys124-Gly137 Purity: >90% Traits: Freeze-dried powder. Buffer formulation: PBS, pH7.4. Applications: Immunogen; Coating Antigen; Blocking Peptide; ELISA.

#### [ <u>USAGE</u> ]

Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

#### [STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.


**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.


### [SEQUENCE]

KDHPLTC DHPQLSG

## Coud-Clone Corp.

## [IDENTIFICATION]








Figure 2. MASS Spectrometry of the Original Peptide

# Coud-Clone Corp.




Figure 3. SDS-PAGE of OVA Conjugated LH