

RPA372Mu01 5mg

Recombinant Gelsolin (GSN)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Ala432~Ser625

Tags: N-terminal His Tag

Subcellular Location: Secreted, Cytoplasm

Purity: > 95%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% Sarcosyl, 5% Trehalose.

Original Concentration: 600µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.0

Predicted Molecular Mass: 24.2kDa

Accurate Molecular Mass: 27kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```

                                AAQHGMDDD GTGQKQIWRI
EGSNKVPVDP ATYGQFYGGD SYIILYNYRH GGRQGQIIYN WQGAQSTQDE
VAASAILTAQ LDEELGGTPV QSRVVQGKEP AHLMSLFGGK PMIIYKGGTS
RDGGQTAPAS IRLFQVRASS SGATRAVEVM PKSGALNSND AFVLKTPSAA
YLWVGAGASE AEKTGAQELL KVLRS
  
```

[IDENTIFICATION]

TGGGCTCAGCAGGCAATGGATGATGGAACCTGGGCGAAGCAGATCTGGAGATGAGGTTCGACAGGTCCAGTGGACCTGCGGATACGGGCGATCTCTATGAGGCGACAGCTACATCTCTGTACACTACGGGCAAGGTGGGCGCTCAGGSHCAGATCTCTACACTGGCAGGCTGCGGAGTCTACCGAGGATGAGGTTC
 A A Q H G H D D G T G Q K Q I W R I E G S N K V P V D F A T Y G Q F Y G G D S Y I I L Y N Y R H G G R Q G Q I I Y N W Q G A Q S T Q D E V A




Figure. Gene Sequencing (Extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.