

RPF150Ra01 50µg

Recombinant Ubiquitin Specific Peptidase 10 (USP10)

Organism Species: *Rattus norvegicus* (Rat)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Met1~Leu794

Tags: N-terminal His Tag

Subcellular Location: Nucleus, Cytoplasm

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: 10mM HEPES, 150mM NaCl, 3mM EDTA, 0.05% Tween20 (pH 7.4), 5% Trehalose.

Original Concentration: 100µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 4.8

Predicted Molecular Mass: 91.0kDa

Accurate Molecular Mass: 91kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MALHNPQYIFGDFSPDEFNQFFVTPRSSVELPPYSGTQCGIQAEELLDGQEHQRIEFGVDE
 VIEPSDGLQRAPSYSISSTLNPQAPEFILGCPTSKKTPDDIEKDETYSSIDQYPASALALESSNA
 EAETLEND SGAGGLGQ RERKKKKRPPGYYSY LKDGSEEGASPAALVNGHATSVGTNSEGVE
 DPEFMVDMLPSVMPRTCDS PQNPMDLISDPVPDSPFPRTLGGDARTAGLPEGCRETDFEQ
 PCLPTDNLLRTAVTQPNAGADTTENLAVANGKILES LGEGTAAANGVELHTDESADLDPAPK
 ESQSPAESALSVSGAISISQPAKSWASLFHDSKPSSSSPVAYVETKCSPPVPSPLASEKQMEVK
 EGLVPVSEDPVAIKIAELLETVTLVHKPVSLQPRGLINKGNWCYINATLQALVACPPMYHLMK
 FIPLYSKVQRCTSTPMIDSFVRLMNEFTNMPVPPKPRQALGDKIVRDIRPGAAFEPTYIYRLL
 TVIKSSLSEKGRQEDAE EYLGFI LNLHEEMLSLKLLSP THEKHSVSNPGSHLIEDEELEDTG
 EGSEDEWEQVGPKNKTSVTRQADFVQTPITGIFGGHIRSVVYQSSKESATLQPFFTLQLDI
 QSDKIRTVQDALESLVARES VQGYTTKTKQEVEVSRRVTLEKLPPVLVLHLKRFVYEK TGGCQ
 KLVKNIEYPVDLEISRELLSPGVKNKNFKCHRTYRLFV VYHHGNSATGGHYTTDV FQIGLNG
 WLRIDDQTVKVINQYQVVRPSADRTAYLLYRRVDLL

[IDENTIFICATION]




Figure . Gene Sequencing (extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was

used in clinical diagnostic or any other procedures.